

Product Name:	Cell Strainer
Product Application:	NEST cell strainer can be used to obtain uniform single cell suspension which is compatible with 50ml centrifuge tubes of most popular brands on the market. • Sample filtration • Sample preparation for Flow Cytometer • Filtration of viscous proteins from inactivated serum • Isolating primary cells to consistently obtain a uniform single cell suspension from tissues • Preparation of single-cell suspensions of blood cells in bone marrow, pancreas, thymus, and lymph nodes • Stem cell preparation • Preparation of samples prior to freezing
Product Materials:	PP frame + Nylon strainer
Product Shelf Life:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good)
Sterilization:	Yes. Sterilized by E-beam, Sterility Assurance Level: SAL=10⁻⁶. The product has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137
Pyrogens:	Non-Pyrogenic
RNase/DNase:	DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Packaging inspection: qualified
Features	1. Fits nearly all 50ml conical tubes. 2. Tab molded within upper ring enables easy sterile handling. 3. Individually Packaged. 4. Non-pyrogenic and DNase/RNase free, non-cytotoxicity.

Head Office

Email: info@cell-nest.com
 Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
 NEST Scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Product Range


Cat. No.	Spec(um)	Color	/Pack	/Case
258366	40	Blue	1	50
258365	70	White	1	50
258364	100	Yellow	1	50
258369	40	Blue	1	50
258368	70	White	1	50
258367	100	Yellow	1	50

Product picture


Head Office Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Technical Drawing



258366 / 258365 / 258364



258369 / 258368 / 258367

Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)