

Product Name:	Cover Glasses
Product Application:	Suitable for a variety of cell staining, including immunohistochemistry and immunofluorescence staining followed by laser confocal observation.
Product Materials:	High quality borosilicate glass
Product Shelf Life:	Temperature range: stored at room temperature Shelf life: 3 Years after date of production (ensure package is in good)
Sterilization:	Yes.
Pyrogens:	Non-pyrogenic
RNase/DNase:	DNase/RNase free
BSE/TSE:	These products are deemed animal free.
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Cell culture test: qualified Packaging inspection: qualified
Features	1. High quality borosilicate glass, resistant to organic solvents. 2. No pyrogen, no cytotoxicity, DNase/RNase free 3. Autoclavable 4. Not TC treated. 5. Thickness: 0.16±0.02mm 6. Can be used in confocal experiment. 7. Each package has an independent batch number for easy quality tracking and traceability.
Packages	Each package has a spare quantity to balance the loss during transportation and use.

Product Range

Cat. No.	Glass Diameter (mm)	Thickness (mm)	Application	/Pack
801010	14	0.16±0.02	24 Well Cell Culture Plate	100
801007	15	0.16±0.02	24 Well Cell Culture Plate	100
801011	18	0.16±0.02	12 Well Cell Culture Plate	100
801008	20	0.16±0.02	12 Well Cell Culture Plate	100
801009	25	0.16±0.02	6 Well Cell Culture Plate	100


Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)

Product Picture



Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Oversea

NEST USA (New Jersey/ Phoenix)
NEST Scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
NEST Scientific (MENA) FZE (Sharjah, United Arab Emirates)