

ISO9001 & ISO13485 & ISO11137

Product Name:	Carboy	
Product Application:	1.Storage container for raw pharmaceutical materials or culture media that require autoclaving sterilization. 2.Storage container for bulk raw pharmaceutical materials or other substances. 3.Storage container for sterile water.	
Product Materials:	PP (barrel/tap) + PP.TPE (cap) <i>meets USP, Class VI standards.</i>	
Product Specifications:	Temperature range: stored at room temperature. Shelf life: 3 Years after date of production	
Sterilization:	NO	
Pyrogens:	Non-Pyrogenic	
RNase/DNase:	DNase/RNase free	
BSE/TSE:	These products are deemed animal free.	
Performance Testing:	Each manufacturing lot is sampled and tested in accordance with standard operating procedures. Appearance inspection: qualified Sealing test: qualified Solvent storage test: qualified Autoclave test: qualified	
Features	1.It is autoclavable. 2.Large capacity solvent storage. 3.Good material resistance and a wide range of solvents can be stored. 4.The sealing performance is secured by the TPE gasket and the thread on the finish, which is well matched to the cap. 5. With the faucet installed, it's easy to drain solvents.	

Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Overseas

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137



Specifications

Volume(L)	With handle	Sterilization	/Case	With tap	Without tap
10	Yes	NO	4	789001	/
10	Yes	NO	6	/	789011
20	Yes	NO	3	789101	789111

Cat. No.	Product Description	/Pack	/Case
789901	83B 3-Port Cap, 1/4"ID Inlet and Outlet Ports, 1/4" Vent Filter Port, Non-Sterile	1	2
789911	83B 3-Port Cap, 1/2"ID Inlet and Outlet Ports, 3/8" Vent Filter Port, Non-Sterile	1	2

Technical Drawing

10 L Carboy, Autoclavable, with Handle, with Tap



Head Office


Email: info@cell-nest.com
 Online: www.cell-nest.com

Overseas


NEST USA (New Jersey/ Phoenix)
 NEST scientific 株式会社 (Yokohama, Japan)
 NEST Scientific Europe B.V (Netherlands)
 Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137

20 L Carboy, Autoclavable, with Handle, with Tap



10 L Carboy, Autoclavable, with Handle, without Tap



Head Office


Email: info@cell-nest.com
Online: www.cell-nest.com

Overseas


NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)

ISO9001 & ISO13485 & ISO11137


20 L Carboy, Autoclavable, with Handle, without Tap



789901



789911



Head Office

Email: info@cell-nest.com
Online: www.cell-nest.com

Overseas

NEST USA (New Jersey/ Phoenix)
NEST scientific 株式会社 (Yokohama, Japan)
NEST Scientific Europe B.V (Netherlands)
Nest Scientific (MENA) FZE (Sharjah, United Arab Emirates)